




Femoral Impactor with Slap Hammer

Providing femoral trial extraction as well as impaction of final components.

The Enztec Femoral Impactor with Slap Hammer is efficient and easy to use, and can be configured to work with any implant system. The integral slap hammer handle allows for efficient removal of femoral trials whilst still allowing impaction with a mallet for implant seating.

With a robust design, you can rely on the Enztec Femoral Impactor with Slap Hammer to support placement, impaction and trial extraction, while ensuring stability and protection of femoral components.



Compatible implants

- Medacta GMK
- United U2
- + more

Features

- Intuitive to use
- Integral slap hammer handle
- Configurable to suit your femoral components
- Easy to clean
- Replaceable impaction head

Specifications

- Materials: Stainless steel, Propylux
- Weight: 910g / 2.0lbs
- Length: 235mm / 9.3in
- Add your logo: available