




# Offset Cup Impactor X

## A new level of in-situ visualization.

With its detachable trinket the Enztec Classic Offset Cup Impactor X can be quickly and easily removed for visual inspection in-situ or intraoperative imaging to verify cup position before final impaction.

The generous offset supports a variety of patient anatomies meaning this instrument can be used in any approach. It is already compatible with a wide range of acetabular cups or will be adapted to fit your cup.

You can rely on the Offset Cup Impactor X to give you the flexibility and confidence to seat your cup in any situation.



### Compatible implants

- Smith + Nephew R3
- DePuy Pinnacle
- DePuy Emphasys
- Stryker Trident
- Signature Logical G
- Corin Trinity
- + more

### Features

- Allows in-situ visualization and intraoperative imaging of the cup before final impaction
- Infinitely adjustable cup orientation for lag screw positioning
- Release lever clear of the incision

### Specifications

- Materials: Stainless Steel, Propylux, Silicone
- Weight: 1300g / 2.9lbs
- Length: 420mm / 16.5in
- Offset: 55mm / 2.2in
- Add your logo: available