





Adjustable Tibial Stylus

Tailored, controlled tibial resection.

With its lightweight all-in-one design, the Enztec Adjustable Tibial Stylus provides flexibility and control to optimize resection depth for your patients.

The rotating and translating stylus tip offers a reliable reference to the tibial plateau, while the easy grip dial allows precise adjustments in resection depth. Whether you require a bone-preserving 1mm adjustment or a more substantial 12mm, you can effortlessly set the resection depth to align with your patient's needs and implant specifications.

The Enztec Adjustable Tibial Stylus gives you tailored resection depth control.



Features

- Customizable resections from 0 to 12mm
- Rotating and translating stylus tip for accuracy
- Easy grip dial for fine adjustments
- Compatible with a 1.37mm saw slot
- Open design for enhanced cleaning and sterilization

Specifications

- Materials: Stainless Steel
- Weight: 160g / 0.4lbs
- Length: 130mm / 5.0in
- Add your logo: available