



Quick Connection Handles

Lightweight, quick connection handles that can be added to any set.

These handles serve as user-friendly tools that simplify the process of establishing connections between devices, ensuring swift and hassle-free interactions.

Choose from the quick release T-Handle or the Torque Indicating handle.



Features

- Comfort and slip resistant silicone handles
- Torque indicating handle indicates up to 4Nm of torque applied
- Light weight and easy to use
- T handle available with AO or Zimmer coupling
- Torque indicating handle available with AO coupling.

Specifications

- Materials: Stainless Steel, Silicone
- Add your logo: not available
- Enquire for more detailed dimensions and weights