




Modular Anvil

Ergonomic handle for modular workflows.

The Enztec Modular Anvil is designed for seamless compatibility across Enztec's modular platform, including the Modular Offset Cup Impactor and the Modular Broach Handle, as well as deSoutter IMT-style connections.

Ergonomically shaped for comfort and secure grip, it is available in Square or Offset configurations to support different workflows and user preferences.

Simple and intuitive to use, the Modular Anvil connects and disconnects quickly and securely, helping maintain procedural efficiency.



Features

- Ergonomically shaped handle
- Recessed, easy to actuate release levers
- Available in Square or Offset
- Orientation agnostic to suit individual workflows

Specifications

- Materials: Stainless Steel
- Weight: 680g / 1.5lbs
- Length: 137mm / 5.4in
- Add your logo: available