



Acetabular Reamers

Reliable cutting performance for accurate repeatable performance.

Enztec Acetabular Reamers are engineered to create a perfectly hemispherical cavity with a smooth finish, supporting precise sizing and confident preparation in every procedure.

Designed for consistent performance, they produce a clean, reliable cut that surgeons can depend on.

Available in 1 mm increments from 36 to 74 mm, with crossbar or bridgeback connections, there is a reamer to suit your needs.





AVAILABLE IN:
- Crossbar
- Bridgeback

Features

- Sharp cutting edge
- Accurate and consistent sizing
- Wide range of sizes to support different patient anatomies (36mm to 74mm)
- Available in 1mm size increments
- Crossbar or bridgeback connection

Specifications

- Materials: Stainless Steel
- Weight: 22-42g (each)
- Diameter: 36mm to 74mm