



# Straight Reamer Driver

Simple, reliable and easy to clean.

Offered in long and short variants, the Enztec Straight Reamer Driver meets the demands of a variety of surgical approaches and patient anatomies. With four power couplings and compatibility with Crossbar and Bridgeback reamers, there is a configuration to suit your needs.

Designed and manufactured with reliability in mind, it features a removable sleeve to protect soft tissue during reaming and to allow easy cleaning. It can be used in forward and reverse reaming, allowing impaction of autologous bone graft for defects.

With a low profile and clean lines, the Enztec Straight Reamer Driver is a solution that supports minimally invasive techniques without compromising performance.



SHORT LENGTH 220MM



STANDARD LENGTH 310MM



## Features

- Removable soft tissue protection sleeve
- Reamer capture allows forward and reverse reaming
- Compatible with Crossbar or Bridgeback reamers
- Standard or Short length
- Easy to clean

## Specifications

- Materials: Stainless Steel, Propylux
- Standard - Weight: 285 - 340g / 0.6 - 0.7lbs  
- Length: 310mm / 12.2in
- Short - Weight: 220 - 250g / 0.5 - 0.6lbs  
- Length: 220mm / 8.7in
- Add your logo: available