


# Adjustable Distal Femoral Alignment Guide

Femoral Resection, aligned to your approach.

The Enztec Adjustable Distal Femoral Alignment Guide offers flexible configurations to support femoral bone cuts aligned to native patient anatomy, including kinematic alignment.

The guide consists of two adjustable paddles, offering the ability to compensate for cartilage that is worn or no longer present, with optional visual feedback on varus/valgus angle through the sagittal alignment connector if desired.

This device is tactile, easy to use, and offers an open design for easy cleaning and sterilisation.



## Features

- Supports femoral bone cuts aligned to native anatomy
- Paddles with 0-4mm range to compensate for cartilage wear
- Intuitive to use
- Modular system configurable to suit your needs
- Full disassembly to optimize tray space
- Open design for easy cleaning.

## Specifications

- Materials: Stainless Steel
- Add your logo: not available
- Enquire for dimensions and weights