

Single Use Pins

Enztec single-use pins, available in smooth, threaded and fluted variants, are the perfect companion for cutting blocks, guides and arrays.

Available in sterile, or non-sterile configurations and made from 455 stainless steel, these pins are high quality ensuring consistent performance without compromise.

Enztec's single use pin range is expanding. If you don't see what you are after please get in touch to discuss your specific pin requirements.


SMOOTH



FLUTED



THREADED



Features

- Tri-flat drive connection
- Trocar-tip
- Sterile or Non-Sterile

Specifications

- Materials: Stainless Steel
- Diameter: 3.2mm
- Length Options: 70mm, 110mm
- Weights:
 - Smooth: 4g (70mm), 6.5g (110mm)
 - Fluted: 3.9g (70mm), 6.2g (110mm)
 - Threaded: 3.9g (70mm), 6.2g (110mm), 8.1g (150mm)