

Angled Reamer Driver




# Angled Reamer Driver

The Angled Reamer Driver supports a Direct Superior Approach with a 55° offset

## Features:

- Connects to industry standard and Enztec crossbar reamers
- Unique capture allows reamers to be securely held but not locked in place
- Disassembles for easy cleaning
- Field changeable drive coupling
- Intuitive assembly
- Replaceable parts, including drive train
- Index-able (45° increments) handle.



Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Angled Reamer Driver, visit:  
[www.enztec.com](http://www.enztec.com)  
[enquiries@enztec.com](mailto:enquiries@enztec.com)  
NZ +64 3 348 0203

**enztec**