


Classic Patella Clamp

Classic Patella Clamp

The Classic Patella Clamp allows quick and efficient patella preparation and cement pressurisation. Its light weight design allows the clamp to be left in place unsupported with minimal risk of patella tendon damage.

Features:

- Light weight
- Fast clamping/releasing
- Single and multisize drill guides available
- Annular rings on multisize drill guides ensure central patella location
- Modular design allows a single handle to be utilised for both drilling and clamping
- Works with both round and anatomic patella.



Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Classic Patella Clamp, visit:
www.enztec.com
enquiries@enztec.com
NZ +64 3 348 0203

Enztec