

Classic Patella Saw Guide




Classic Patella Saw Guide

Where simplicity meets functionality – this intuitive ratcheting mechanism ensures secure and accurate patella resection.

Features:

- Accommodates different resection heights (8, 9, 10mm) through a fixed three-level stylus
- Single-handed operation to engage and lock onto the patella
- Simple construction – no removable parts
- Easy access for cleaning
- For use with a 1.27mm saw blade.



For further information or to order your Enztec Classic Patella Saw Guide, visit:

www.enztec.com

Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

enquiries@enztec.com | NZ +64 3 348 0203

 **enztec**