


Charnley Retractor Suite

A self-retaining frame and retractor system utilized during total hip arthroplasty, enabling a clear view and retraction of the hip capsule and musculature.

Enztec Charnley Standard Frame and Retractors

PRODUCT DESCRIPTION	PRODUCT CODE	A LENGTH	B LENGTH
Standard Frame	4050-0150	256mm	200mm
Clawed Retractors:	50mm 4050-0161	170mm	50mm
	75mm 4050-0162	170mm	75mm
	100mm 4050-0163	170mm	100mm
	125mm 4050-0164	170mm	125mm
Spiked Retractors:	50mm 4050-0171	170mm	50mm
	75mm 4050-0172	170mm	75mm
	100mm 4050-0173	170mm	100mm
	125mm 4050-0174	170mm	125mm

*Enztec Charnley Retractors are compatible with both standard and fixed frames.



For further information or to order your Charnley Retractor Suite, visit: www.enztec.com

Enztec is an MDR-certified, ISO 13485 designer and manufacturer of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries. enquiries@enztec.com | NZ +64 3 348 0203




Charnley Retractor Suite

A self-retaining frame and retractor system utilized during total hip arthroplasty, enabling a clear view and retraction of the hip capsule and musculature.

Enztec Charnley Fixed Blade Frame and Retractors

PRODUCT DESCRIPTION	PRODUCT CODE	A LENGTH	B LENGTH
Fixed Blade Frame	4050-0151	256mm	200mm
Clawed Fixed Blade Retractors:	50mm	4050-0181	79mm
	75mm	4050-0182	79mm
	100mm	4050-0183	100mm
	125mm	4050-0184	79mm
Spiked Fixed Blade Retractors:	50mm	4050-0191	79mm
	75mm	4050-0192	79mm
	100mm	4050-0193	100mm
	125mm	4050-0194	79mm



For further information or to order your Charnley Retractor Suite, visit: www.enztec.com

Enztec is an MDR-certified, ISO 13485 designer and manufacturer of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries. enquiries@enztec.com | NZ +64 3 348 0203

