



Femoral Impactor with Slap Hammer



# Femoral Impactor with Slap Hammer

---

The Femoral Impactor with Slap Hammer allows for femoral trial extraction in addition to femoral component insertion, building on the proven functionality of our ever popular femoral impactor and making it a valuable addition to any primary or revision knee set.

## Features:

- Simple Slap Hammer that telescopes out for femoral trial extraction
- Screwable impaction handle allows a wide range of femoral components to be firmly gripped
- A range of standard arm tips available to ensure a secure hold of femoral components
- Impactor head is replaceable if required due to cement deposit, adhesion etc
- Contoured arms for easy grip.



Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Femoral Impactor with Slap Hammer, visit: [www.enztec.com](http://www.enztec.com)  
[enquiries@enztec.com](mailto:enquiries@enztec.com)  
NZ +64 3 348 0203

**Enztec**