




# Classic Patella Clamp

For quick and efficient patella preparation and cement pressurization.

Light-weight and modular, the Enztec Classic Patella Clamp quickly transforms from drilling to clamping modes. Designed with ease-of-use in mind, its light weight means it can be left unsupported during cement pressurisation with no risk of damage to the patella tendon.

The Enztec Classic Patella Clamp is versatile and can be configured to suit your needs.



## Features

- Modular attachments
- Single and multisize drill guides tailored to your implant
- Fast clamping and releasing
- Can be left in place unsupported with no risk of patella tendon damage
- Compatible with anatomic and round patella

## Specifications

- Materials: Stainless Steel, Propylux, Silicone
- Handle weight: 315g / 0.7lbs
- Handle length: 220mm / 8.7in