




# Femoral Head Impactor

A reliable tool for controlled impaction.

With a deep dish the Enztec Femoral Head Impactor ensures stable location and elimination of slippage during impaction.

Designed with a short length, this instrument gives great impaction control. With a removable quick connect head, it is also easy to clean.

You can rely on the Enztec Femoral Head Impactor to give you stable, controlled impaction.



## Features

- Easy to clean removable head that can be quickly reconnected
- Compatible with 28, 32, 36 and 44mm femoral heads
- Replaceable head
- Short length for great impaction control

## Specifications

- Materials: Stainless Steel, Propylux
- Weight: 390g / 0.9lbs
- Length: 160mm / 6.3in
- Add your logo: available