




Premium Offset Cup Impactor

Quick cup release for maximum intra-operative efficiency.

The Enztec Premium OCI facilitates easy lag screw positioning and quick withdrawal from the cup. With a minimal offset design, the Enztec Premium OCI is optimised for MIH procedures and can be used in any approach.

Clinically proven with over 10 years use in the field, the Enztec Premium OCI is your go to tool for efficient cup placement and impaction.

The Enztec Premium OCI is compatible with a range of acetabular cups.



Compatible implants

- B-one
- Amplitude C2
- Smith + Nephew R3

Features

- Easy lag screw positioning
- Quick release mechanism for efficient detachment from the cup in situ
- Full disassembly for cleaning

Specifications

- Materials: Stainless Steel, Propylux
- Weight: 1100g / 2.4lbs
- Length: 420mm / 16.5in
- Offset: 25mm / 1.0in
- Add your logo: available