



Premium Pin Puller

A high end solution for discerning surgeons offering quick, fluid extraction of pins with minimal effort.


Compressing the handle grips and extracts the pin in one fluid motion. The jacking motion results in efficient, effortless pin removal, reducing operator strain and effort.

Simple and easy to use, the Enztec Premium Pin Puller gives you confidence to extract pins in the most difficult circumstances.



01. LOCATE OVER PIN

Low profile nose aids ease of use



02. COMPRESS HANDLE

Compressing the handle grips and extracts the pin in one easy motion



Features

- Compatible with 3.2mm / 0.125in Steinmann pins
- Effortless pin extraction
- Contoured handle provides comfort and grip

Specifications

- Materials: Stainless Steel
- Weight: 220g / 0.5lbs
- Length: 155mm / 6.1in
- Add your logo: available