


# Straight Reamer Driver

Simple, reliable and easy to clean.

With four different power couplings available and compatible with Crossbar or Bridgeback reamers, there is an Enztec Straight Reamer Driver to suit your needs.

Designed and manufactured with reliability in mind, the Enztec Straight Reamer Driver features a removable sleeve to protect soft tissue during reaming and allow easy cleaning. It can be used in forward and reverse reaming, allowing impaction of autologous bone graft for defects.

With a low profile, the Enztec Straight Reamer Driver supports minimally invasive techniques.



## Features

- Removable soft tissue protection sleeve
- Reamer capture allows forward and reverse reaming
- Compatible with Crossbar or Bridgeback reamers
- Fits in most instrument trays (standard length and profile)
- Easy to clean

## Specifications

- Materials: Stainless Steel, Propylux
- Weight: 285 - 340g / 0.6 - 0.7lbs
- Length: 310mm / 12.2in
- Add your logo: available