

PCL Retractor




# PCL Retractor

---

The Enztec PCL Retractor sits in the femoral condylar notch and straddles the cruciate ligament, providing better access to the joint.

Efficient and comfortable to use, Enztec retractors minimise soft tissue damage and are readily available at a competitive price.



For further information or to order your Enztec PCL Retractor, visit:  
[www.enztec.com](http://www.enztec.com)

Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.  
[enquiries@enztec.com](mailto:enquiries@enztec.com) | NZ +64 3 348 0203

