


Premium Patella Saw Guide

Premium Patella Saw Guide

The Premium Patella Saw Guide securely grasps the patella on the medial-lateral boundaries through an ergonomic, single-handed operation. This allows accurate patella resection, ensuring correct flexion and extension of the leg can be obtained.

Features:

- Two variants available offering a 9mm or 10mm standard resection
- Bone preservation allowed for through the use of optional 1mm, 2mm or 3mm shims
- Fine tooth profile with mixture of long/short teeth provides positive engagement with patella
- Pistol grip for ergonomic handling
- Dual release lock buttons prevent accidental release
- Finger grip located on jaw allows for forceful release, should teeth stick in patella
- For use with a 1.27mm saw blade.



Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Premium Patella Saw Guide, visit:

www.enztec.com
enquiries@enztec.com
NZ +64 3 348 0203

The logo consists of a stylized 'E' symbol followed by the word 'enztec' in a bold, lowercase sans-serif font.