



Classic Pin Puller

Classic Pin Puller


The Classic Pin Puller is a sleek instrument designed to compliment your instrument sets. Securely gripping pins with minimal length exposed, this instrument is simple and easy to maintain.

Features:

- For removal of 1/8" (3.2mm) Steinman Pins
- Ergonomic handle design
- Can extract pins of any length
- Able to grip a pin with only 1/8" (3.2mm) exposed
- Patented design (US15/004482, USD796,039)

Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Classic Pin Puller, visit:
www.enztec.com
enquiries@enztec.com
NZ +64 3 348 0203



enztec