

Classic Pin Puller with Anvil

Classic Pin Puller with Anvil

The Classic Pin Puller with Anvil provides additional support for a comfortable, non-slip grip, while also providing an impaction face if required.


Features:

- Able to grip a pin with only 1/8" (3.2mm) exposed
- For use with 1/8" (3.2mm) Steinman Pins
- Ergonomic handle design
- Can extract pin of any length.

Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your Enztec Classic Pin Puller with Anvil, visit:

www.enztec.com
enquiries@enztec.com
NZ +64 3 348 0203



ENZTEC