




Offset Reamer Driver – Open

# Offset Reamer Driver – Open

The Open Offset Reamer Driver allows the reamer to be placed in the incision prior to connection, supporting a Direct Anterior Approach.

## Features:

- Can be detached from the reamer in situ, allowing quick visualization of the ream
- Connects to industry standard and Enztec crossbar reamers
- Unique capture mechanism allows reamers to be securely held but not locked in place
- Index-able handle (45° increments)
- Field changeable drive coupling
- Disassembles for easy cleaning
- Replaceable parts, including drive train
- Fits in most existing reamer pans with industry standard length and profile
- Intuitive assembly.



Enztec is focused on the design and manufacture of orthopaedic instruments. We are a proven, quality approved supplier to all major orthopaedic implant companies and welcome your enquiries.

For further information or to order your **Enztec Offset Reamer Driver – Open**, visit: [www.enztec.com](http://www.enztec.com)  
[enquiries@enztec.com](mailto:enquiries@enztec.com)  
NZ +64 3 348 0203

**Enztec**